

Testing

• How Participants Test

Someone is considered eligible for testing if they are: symptomatic, exposed, or high risk immunocompromised, or pregnant.

STEP 1: Registration

- Participants register their test using the QR codes or URLs on the 11x17 signage or registration fliers. Please hand out a registration flier along with each test.

STEP 2: Activating a Color Kit

- Participants must activate a sample in order for it to be processed. If a sample arrives at the lab unactivated, it will be destroyed and will not be processed.

Video:

[Registration and sample activation](#) →

Enter your barcode information.

Locate the barcode (starts with the letter D) and accession number (starts with the letter C) on your collection tube. Enter them below. Having trouble? [Get help.](#)

Barcode:

Accession number:

[Take or Upload Photo](#)

[Continue](#)



Scan for video

STEP 3: Sample Collection

- Participants will collect their sample by following the swabbing instructions in their test kits and on the site set up signage.
- Participants place the tube back into the bag and seal it
- Participants return the sample to the dedicated collection bin
- Participants should take a picture of their barcode card located in the bag with their test kit in case they need to get in touch with Patient Support.

More info:

[General onsite sample collection](#) →



Scan for video

Participants must complete each step, register, activate and return their sample on the same day, or their test will not be processed and they will not receive results.

Test Results

- **Notifications**

Participants will be notified via text and email (if provided) when their results are ready to be viewed. The text message will not contain information about the type of result (positive, negative or inconclusive).

STEP 1: Results Are Ready Notifications

Participants will receive a text message and email (if provided email address during registration) notifying them that their results are ready.

STEP 2: View Results

Participants click on the notification from the text or email link. The participant will arrive at a home screen that indicates the results are ready and what date the testing was conducted. The most recent test will appear on this home screen.

Participants will need to enter their date of birth and barcode card to access their results.

Participants should take a picture of their barcode card located in the bag with their test kit, in case they need to get in touch with Patient Support. If a participant loses their barcode, please contact Color Support for assistance at 844-352-6567 or mycovidtest@color.com.

